


# High-Density VEP Set-Up Checklist

## MONITOR

- Turn on the **CRT monitor** on the cart in 120G

## AMPLIFIER

- In 120G Chandlee remove cables from amplifier and connect them to pegs on wall.



- Wheel amp between the monitor and the subject chair so that it is to the right of the subject chair. **ALWAYS PULL AMP BACKWARDS TO PREVENT TIPPING**
- Plug the 4 cords from the wall ports connecting rooms G and H into respective ports (Digital input, USB, and TTL input external timebase) of the amplifier.
- Turn power on (I = on).



- Turn small switch on amplifier to **NORMAL**
  - If this is on another setting during calibration and impedance check, you will see an error stating that everything is out of specification.
  - If there is an error, you must begin the setup from step one again.

## COMPUTERS

- In 120H Switch on power of large surge protector on bottom left shelf.

- Switch on power of smaller, thin surge protector on second shelf (ignore lights; always on).
- Turn on **PD Video** and **NetStation** computers, and both monitors.
  - If either monitor screen is unresponsive, you may need to toggle to the correct computer via the gray KVM switches on the top left (PD Video / Link 15) and right (NetStation / PD Host) shelf.
- Turn on external timebase, behind the monitors (switch is back, right side).
- Remove the plug from the external timebase outlet labeled **OUT**. (Clock Out stays in).




- Make sure that white round switch on top of the left KVM switch is set to B.



- On PD Video computer, open Power Diva Video 3.4 software. (icon label says alias)



- In Power Diva Video, go to Configuration > Video Manager. In **Mode** window, make sure the calibrated video setting is selected:
  - 800 x 600, 72 Hz, 8 bit



- Click OK

## Prepare Stimuli

### For Child and Adult

- In Run mode (pull-down tab), choose **Auto Sequence**.
- In Msg mode (pull-down tab), choose **Net Station**.
- In Stim Set, Go to Open, and choose the appropriate stimulus set.
  - Child
    - mofo\_tuning\_child\_20150428
  - Adult
    - moco\_3pat\_hidens\_20150428



### For Infant

- In Run mode (pull-down tab), choose **Baby Mode**.
- In Msg mode (pull-down tab), choose **Net Station**.
- In Stim Set, Go to Open, and choose the appropriate stimulus set.
  - inf\_2pat\_lamrad\_20150428 (most used)
  - inf\_2pat\_lamrot\_20150428
  - inf\_2pat\_radrot\_20150428

Stim Set: inf2patlamrad [New] [Open] [Close]

Stimulus: 1 MOCO stim [New] [Clone] [Edit] [X]

Run Mode: Baby Mode [Run] [Reset]

Msg Mode: NetStation ☐ Log Response

Status: Ready

- Press the **Run** button. In the DAQ Rate Settings window, there is a binary code to confirm. On the side of the external time box facing you, there are silver switches (**up = 1, down = 0**). Read the binary code in the PD Video window from **RIGHT to LEFT**, and check to make sure that the switches on the timing box correspond to this code. If they don't match, notify the PI. If they do, press **Cancel**.

### AutoRun Settings

**Response Log File**

Subject Name:

Log File Name:

You have requested recording of the subject's responses. Please provide the Subject's name. It will be used as part of the name of the Log File as shown above.

**DAQ Rate Settings**

DAQ Rate (Hz): 432.43 [Binary Code: 000001101111]

Choose the frequency to be used as the external clock rate for NetStation data acquisition, and dial corresponding binary code on the switch box.

**NetStation Recording Control**

☐ Start/Continue NetStation recording

By checking this box you request turning on the NetStation recording mode every time the stimulus presentation begins.

To create a new recording use the NetStation's "Waveform Recorder Controls" panel.

**Run Time Segment**

Concurrent Task Mode: No Task

Resp Dev: Keyboard

Conditions: 1 MOCO Target Param:  Set vs

Staircase: 2D1U 1-2, 82% Letter Scale: 1 N Symb: 1

Probe Duration min/max (sec):   Set Task

[Cancel] [Run]

## Prepare Net Station

- Go to the Net Station computer and log in with the **gilmorelab** account (w/ lab password - **0baby0**).

- In the dock at the bottom of the screen, click on the icon for Net Station (purple icon) to open the program.



- When the window pops up, there should be a key symbol with the words **Wenger** on the bottom left-hand side. If it is not there, this means that the NS HASP key is not plugged in.
- In the Net Station window, click **New Session**
- Select the session template labeled **Power Diva Setup**.

## AFTER PARTICIPANT ARRIVES

- Enter participant ID, Experiment Name, Date of birth, and net serial number. Lab uses the following standard ID number for each session: YYMMDD\_ID# (Test Date\_Participant ID number).

**Enter Session Information**

Identify Subject: Enter Subject Identifier or Select From Table.

20121106\_1200

- 20120706\_0900
- 20120710\_0300
- 20120721\_1100
- 20120724\_1000
- 20120724\_1100
- 20121106\_1200

Session Name: 20121106\_1200 20121106 13.07 Rename...

Session Information:

Experiment ID	moco3pat_hidens		
Date of Birth	08	06	2006
	MM	DD	YYYY
Geodesic Sensor Net Serial #	H1861		

Cancel Session Begin Session


- **Begin Session.**
- Run Zeros and Gains (should run automatically). The net cannot be plugged into the articulating arm at this time.



- Complete [EEG Net Placement](#). (separate protocol)
- After netting, plug the net into the articulating arm.




- Measure Impedance
  - Go to the **Panels** menu at the top of the screen and select **Impedance**
  - Then click the **Measure Impedance** button.
  - Adjust channels as necessary.
  - Click **Save & Close** when satisfactory.
  - If the incorrect impedance check appears, double check that both connectors on the amplifier are securely plugged in.




- Go to the **Panels** menu at the top of the screen and select **Dense Waveform Display**.





- Turn small switch on amplifier to **EXTERNAL TIMEBASE**.



## Experimental Setup

- The monitor should be located 60cm from the bridge of the nose on the participant
- The chair height should be set so the participant is looking directly at the **X** in the middle of the screen.
- High Density
  - Baby Mode (Researcher in the room)
    - Record: Use the Keyboard to record (**R**)
      - When recording there is an **X** in the middle of the screen.
    - Pause Study: Use the keyboard to pause (**W**)
      - When paused there is an **X** with a **circle** around it in the middle of the screen
    - Use the toys to keep attention toward the monitor
  - Auto Sequence
    - Pause Study: The keyboard **space bar** in the testing chamber or on the PD Video computer
    - Continue Study: The **Continue** button must be clicked in the Power Diva software
- Low Density
  - Game controller is used to record and pause

## Run Experiment

- In Power Diva Video, press the **Run** button in the **Stim Set** window
- The **Auto Run Settings** window will open
- **Check the box** in the Net Station Recording control window labeled **Start/Continue Net Station Recording**
- Press **Run**.





Subject Name

Log File Name

You have requested recording of the subject name. It will be used as part of the name of

### DAQ Rate Settings

DAQ Rate (Hz)

Choose the frequency to be used as the external clock and dial corresponding binary code on the s

### NetStation Recording Control


☒ Start/Continue NetStation recording

By checking this box you request turning on stimulus presentation begins.

To create a new recording use the "New Recording" panel.



- Try to get the following number of blocks of data:
  - Infant: 6
  - Child: 10
  - Adult: 12
- On the Net Station computer, make sure the "Recording" symbol is displayed on the upper left-hand side of the screen. Also, make sure the event flag labeled "DIN4" is present, as well as a flag below this with the condition number (e.g., "c004").
  - the flag "p003" means that the study was paused





---

## AFTER SESSION

- Once recording is complete, select the **close session** button on the Net Station computer (this saves the session).
- Return computers and other devices to their default settings.
- Return and replugin amplifier to its default setup.
- Clean net and begin [data export](#) (separate protocols).